



Certified Orthopedics, Inc.
2415 East Mulberry Street, Unit 7
Fort Collins, CO 80524
ph 970-482-7116 • fax 970-498-9529
www.certifiedortho.com

CERTIFIED

DYNA-LOCK NIGHT AFO ORDER FORM

P
A
T
I
E
N
T

First Name: _____
Last Name: _____
Age: _____ Sex: _____ Height: _____ Weight: _____
Diagnosis: _____
Phone #: _____ Fax #: _____
☐ Right ☐ Left ☐ Bilateral

P
R
A
C
T
I
T
I
O
N
E
R

Date: _____ P.O.#: _____
Facility: _____
Address: _____
Contact Name: _____
Ship via: _____
Requested Date of Delivery: _____
(Please allow for shipping time)

S
T
Y
L
E
S

☐ Cast ☐ Scan ☐ Custom to Measurement



☐
DYA-ST1
Custom Dyna-Lock
Night Stretch AFO
Free Motion Medial Joint
☐
DYA-ST2
Custom Dyna-Lock
Night Stretch AFO
w/R.O.M Medial Joint

MODIFYING INSTRUCTIONS

ANKLE ☐ As Is ☐ Correct to 90°
☐ Amount _____° ☐ Dorsi ☐ Plantar
HINDFOOT ☐ As Is ☐ Correct to Neutral
FOREFOOT ☐ As Is ☐ Correct to Neutral
☐ Metatarsal Pad: Amount _____" ☐ S.T. Modification
☐ Other Modifying Instructions: _____
☐ Transfer Pattern #: _____





STANDARD DYNA-LOCK NIGHT AFO INCLUDES:

- Standard Modifications
- 3/16" PolyPro Plastic
- Dyna-Lock Ankle Joint with Medial Free or R.O.M. Joint
- O-Ring Engagement
- Flexible Plastic SMO Lined with Aliplast
- Flexible Plastic Calf Lacer Lined with Aliplast
- Dorsum Strap & Ankle Strap
- Dorsum Overlaps
- Malleolus Pads
- Non-Skid Foot Plate
- Posterior Stop Strap
- Padded Tibial Strap
- Transfer Pattern

EXTERNAL POST

☐ Crepe	☐ Plastic	☐ Other _____
FOREFOOT	HINDFOOT	BOTH
☐ Neutral	☐ Neutral	☐ Neutral
☐ Medial	☐ Medial	☐ Medial
☐ Lateral	☐ Lateral	☐ Lateral
Amount: _____	Amount: _____	Amount: _____

*Non-Loaded & Neutral
*MM Measurements



A. _____ Supramall. Circ.	K. _____ Width Apex of Nav. & Base of 5th Meta
B. _____ Instep Circ.	L. _____ Width of Heel at Widest
C. _____ Midfoot Circ. through Arch	M. _____ Length Heel to Medial Mall
D. _____ 1st Met to 5th Met Circ.	N. _____ Length Ground to Apex of Med. Mall.
E. _____ Foot Length	O. _____ Length Heel to Lateral Mall.
F. _____ Heel to 1st Met Head Length	P. _____ Length Ground to Apex of Lat. Mall.
G. _____ Heel to Navicular Length	Q. _____ Width at Supramalleolar
H. _____ Heel to Base of 5th Met Length	R. _____ Width at Malleoli
I. _____ Heel to Base of 5th Meta. Length	S. _____ Calf Circum. at Widest
J. _____ Width between Met. Heads	T. _____ Height to Fibular Head

ADDITIONAL INFORMATION: _____