



Certified Orthopedics, Inc.
2415 E. Mulberry St., Unit 7
Fort Collins, CO 80524
ph 970-482-7116 • fax 970-498-9529
www.certifiedortho.com

CERTIFIED

STABILIZER LEATHER AFO SYSTEM ORDER FORM

P
A
T
I
E
N
T

First Name: _____
Last Name: _____
Age: _____ Sex: _____ Height: _____ Weight: _____
Diagnosis: _____
Phone: _____ Fax: _____
☐ Right ☐ Left ☐ Bilateral
Cast Scan

P
R
A
C
T
I
T
I
O
N
E
R

Date: _____ P.O.#: _____
Facility: _____
Address: _____
Contact Name: _____
Ship via: _____
Requested Date of Delivery: _____

(Please allow for shipping time)

S
T
Y
L
E
S



FABRICATION INSTRUCTIONS

CLOSURE

- ☐ Laces (standard)
- ☐ Velcro
- ☐ Combination
- ☐ Boot Hooks on top three holes

COLOR

- ☐ Natural Russet
- ☐ Black
- ☐ Brown
- ☐ Taupe


ADDITIONAL OPTIONS

- ☐ Heel Cut Out
- ☐ PTB Height
- ☐ Plastic Reinforced Tongue
- ☐ Sew in Tongue ☐ Medial ☐ Lateral
- ☐ Medial/Lateral Arch Pull Strap
- ☐ Line Inside Hydrostatic with 1/8" Plastazote foam
- ☐ Line Inside Hydrostatic with X-Static
- ☐ Extended Foot Plate - Full Foot
- ☐ Partial Foot Toe Filler on Full Foot Plate
- ☐ Tri-lam Insert ☐ Plastazote ☐ Microcell
- ☐ Extra Padding
- Specify Area: _____
- ☐ Metatarsal Pad
- ☐ Other: _____
- ☐ Correction Flare Padded/Lined
 - ☐ Valgus/Medial ☐ Varus/Lateral

CAST PREPARATION INSTRUCTIONS

- ANKLE** ☐ As Is ☐ Correct to 90°
☐ Amount _____° ☐ Dorsi ☐ Plantar

- HINDFOOT** ☐ As Is ☐ Correct to Neutral
FOREFOOT ☐ As Is ☐ Correct to Neutral



STABILIZER LEATHER AFO's INCLUDE

- Standard Modifications
- Molded leather Cream Cow inside lining
- Plastic reinforced ankle padding with 1/8" Aliplast
- Natural Russet leather lacer
- Lace Closure
- Length of foot plate trimmed proximal to metatarsal heads

ADDITIONAL INFORMATION