



Certified Orthopedics, Inc.
2415 E. Mulberry St., Unit 7
Fort Collins, CO 80524
ph 970-482-7116 • fax 970-498-9529
www.certifiedortho.com

CERTIFIED

STABILIZER DYNAMIC AFO ORDER FORM

P
A
T
I
E
N
T

First Name: _____
Last Name: _____
Age: _____ Sex: _____ Height: _____ Weight: _____
Diagnosis: _____
Phone #: _____ Fax #: _____
☐ Right ☐ Left ☐ Bilateral
Cast Scan

P
R
A
C
T
I
T
I
O
N
E
R

Date: _____ P.O.#: _____
Facility: _____
Address: _____
Contact Name: _____
Ship via: _____
Requested Date of Delivery: _____

(Please allow for shipping time)

S
T
Y
L
E
S



SDA-100
Stabilizer
Dynamic AFO



SDA-200
Dynamic AFO w/
Foot Orthotic
UCBL



SDA-300
Dynamic AFO w/
Molded Inner Boot
LDPE Boot
Flexible Boot



SDA-400
Dynamic AFO w/
Inner Boot & Calf Lacer
LDPE Boot Flexible Boot
LDPE Lacer Flexible Lacer



SDA-500
Dynamic AFO w/
Partial Foot
Transmet
First Ray
Lateral Column

AFO Plastic Choice:

Polypro

Procomp Carbon
Plastic


CAST PREPARATION INSTRUCTIONS

ANKLE ☐ As Is ☐ Correct to 90°
☐ Amount _____° ☐ Dorsi ☐ Plantar
HINDFOOT ☐ As Is ☐ Correct to Neutral
FOREFOOT ☐ As Is ☐ Correct to Neutral
Metatarsal Pad: Amount _____" S.T. Modification
Other Modifying Instructions: _____

ADDITIONAL OPTIONS

Open Heel Trim Line Closed Heel Trim Line
Extended Trim Line to Control Forefoot Medial Lateral
Correction Flare padded/lined Valgus/Medial Varus/Lateral
Proximal Trim Line Flare
Footplate Length Full Foot (Std) Sulcus Proximal Metheads
Add Pretibial Shell
Add Molded Calf Lacer (SDA-100, SDA-200, SDA-300, SDA-500)
LDPE Calf Lacer Flexible Plastic Calf Lacer
Liner: _____ (Specify Area & Material)
Pads: _____ (Specify Area & Material)
Transfer Pattern (specify) _____
Plastic Color (black=standard) _____

ADDITIONAL INFORMATION:



Foot Length: _____

STANDARD DYNAMIC AFO INCLUDES:

- Standard Modifications
- Polypro or Procomp Carbon Reinforced Plastic
- GO AFO Dynamic Carbon Strut
- Ankle Reinforcement
- Custom Foot Orthotic, UCBL, Inner Boot, Calf Lacer, and/or Partial Foot (Depending on Model Selected)
- Desired Trim Lines, Desired Flares, Desired Footplate Length
- Pretibial Strap, Ankle Strap, and/or Dorsum Strap
- Choice of Liner
- Dorsum Overlaps and Malleolus Pads
- Transfer Pattern